



Anti-MSH-alpha polyclonal antibody (DPAB-DC4548)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody may also react with Mouse, Rat and Rabbit
Immunogen	A synthetic peptide (conjugated with BSA) corresponding to MSH-alpha.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr,
Format	Liquid
Size	50 µl
Buffer	In Tris-HCl buffer (1% BSA, 0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	α-Melanocyte-stimulating hormone (alpha-MSH; α-MSH), also known as alpha-melanotropin, alpha-melanocortin, or alpha-intermedin, is a naturally occurring endogenous peptide hormone of the melanocortin family, with a tridecapeptide structure and the amino a
---------------------	---

Keywords

α -MSH; MSH-alpha
