



Anti-MTRNR2L1 (C-terminal) polyclonal antibody (DPAB-DC4549)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids at C-terminus of human MTRNR2L1. The sequence is C-TSEIDL PVKRRT
Source/Host	Goat
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Concentration	0.5 mg/mL
Size	100 µg
Buffer	In 0.5 mg/mL in Tris saline, pH7.3 (0.5% BSA, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	MTRNR2L1
----------	----------

GENE INFORMATION

Entrez Gene ID [100462977](#)

UniProt ID [P0CJ68](#)
