



Anti-Mu Phage Mu Phage gp44 (full length) polyclonal antibody (DPAB-DC4550)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to full length Mu Phage gp44.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Virus
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	25 µg
Buffer	In PBS (0.1% proclin, 2% Block Ace)
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Mu Phage gp44
----------	---------------