



Anti-N. gonorrhoeae Porin Polyclonal antibody (DPAB-DC4560)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	N. gonorrhoeae Porin
Immunogen	Porin from Neisseria gonorrhoeae.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N. gonorrhoeae
Conjugate	Unconjugated
Applications	WB, IF, ELISA,
Format	Liquid
Size	1 mg
Buffer	In PBS, pH 7.4 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Neisseria gonorrhoeae, also known as gonococci (plural), or gonococcus (singular), is a species of Gram-negative coffee bean-shaped diplococci bacteria responsible for the sexually transmitted infection gonorrhea.
--------------	---

Keywords

Porin; Neisseria gonorrhoeae
