



Anti-N. meningitidis type A, B, C Polyclonal antibody (DPAB-DC4562)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	All antigens. Has not been absorbed and may react with related microorganisms. Reactive with serogroups A, B and C.
Target	N. meningitidis type A, B, C
Immunogen	Native from Neisseria meningitidis types A, B and C.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N. meningitidis
Conjugate	Unconjugated
Applications	IF, ELISA,
Format	Liquid
Size	1 ml
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction

Neisseria meningitidis, often referred to as meningococcus, is a bacterium that can cause meningitis and other forms of meningococcal disease such as meningococemia, a life-threatening sepsis. N. meningitidis is a major cause of illness and death during

Keywords

N. meningitidis; Neisseria meningitidis
