



Anti-Nephronophthisis polyclonal antibody (DPAB-DC4564)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant GST fusion protein corresponding to amino acid 612-691 of mouse pic.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB, IP,
Format	Liquid
Size	250 µg
Buffer	In PBS (0.08% sodium azide)
Preservative	0.08% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Nephronophthisis is a genetic disorder of the kidneys which affects children. It is classified as a medullary cystic kidney disease. The disorder is inherited in an autosomal recessive fashion and, although rare, is the most common genetic cause of childh
Keywords	Nephronophthisis (pic)