



Anti-Neurofibrillary Tangle polyclonal antibody (DPAB-DC4567)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	human paired helical filaments from neurofibrillary tangles. Reacts with both intracellular and extracellular neurofibrillary tangles.
Immunogen	Disaggregated paired helical filaments.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr, ELISA, EIA, ID,
Format	Liquid
Size	200 µl
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Neurofibrillary Tangles (NFTs) are aggregates of hyperphosphorylated tau protein that are most commonly known as a primary marker of Alzheimer's Disease. Their presence is also found in
---------------------	---

numerous other diseases known as tauopathies. Little is known about t

Keywords

NFT; Neurofibrillary Tangle
