



## Anti-Neuroketals polyclonal antibody (DPAB-DC4568)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Specificity	Neuroketal adducted proteins.
Immunogen	Neuroketal conjugated.
Source/Host	Goat
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, IHC-P, ELISA, EIA,
Format	Liquid
Size	1 ml
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

### BACKGROUND

Keywords	Neuroketals
----------	-------------