



# Rabbit anti-NFKB1 polyclonal antibody (DPAB-DC4569)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to NFKB1. The sequence is QLRRKSpDLETSEP
Isotype	IgG
Source/Host	Rabbit
Conjugate	Unconjugated
Applications	ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 10 mM PBS
Preservative	None
Storage	Store at -20°C on dry atmosphere. After reconstitution with deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	NFKB1
----------	-------

## GENE INFORMATION

Entrez Gene ID [4790](#)

---

UniProt ID [P19838](#)

---