



Anti-NLRP1 (aa 2-25) polyclonal antibody (DPAB-DC4572)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes NALP1.
Immunogen	A synthetic peptide corresponding to amino acids 2-25 of NLRP1. The sequence is AGGAWGRLACYLEFLKKEELKEFQ
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In buffer containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	NLRP1
----------	-------

GENE INFORMATION

Entrez Gene ID [22861](#)

UniProt ID [E9PE50](#)
