



# Anti-NMU-8 polyclonal antibody (DPAB-DC4573)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Synthetic pig neuromedin U-8 conjugated to BSA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Chicken, Pig, Rat
Conjugate	Unconjugated
Applications	IHC-Fr,
Format	Lyophilized
Size	50 µl
Buffer	Lyophilized from serum
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50-100 uL distilled water, store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	NMU-8
----------	-------

## GENE INFORMATION

Entrez Gene ID [10874](#)

---

UniProt ID [P48645](#)

---