



# Anti-NO-L.Methionine polyclonal antibody (DPAB-DC4582)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                                                                                                                                                                       |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen          | NO-L. Methionine conjugated with carrier protein.                                                                                                                                     |
| Source/Host        | Rat                                                                                                                                                                                   |
| Species Reactivity | N/A                                                                                                                                                                                   |
| Purification       | Ammonium sulfate precipitation                                                                                                                                                        |
| Conjugate          | Unconjugated                                                                                                                                                                          |
| Applications       | IHC,                                                                                                                                                                                  |
| Format             | Lyophilized                                                                                                                                                                           |
| Size               | 100 µl                                                                                                                                                                                |
| Buffer             | Lyophilized from PBS, pH 7.4                                                                                                                                                          |
| Preservative       | None                                                                                                                                                                                  |
| Storage            | Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing. |

## BACKGROUND

|          |                 |
|----------|-----------------|
| Keywords | NO-L.Methionine |
|----------|-----------------|