



Anti-IAV H5N1 Nonstructural Protein 2

Polyclonal antibody (DPAB-DC4588)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Nonstructural protein 2 (NS2) from an H5N1 strain of avian influenza A.
Target	IAV H5N1 Nonstructural Protein 2
Immunogen	A synthetic peptide corresponding to 14 amino acids in internal region region of Influenza A H5N1 (Avian) Nonstructural Protein 2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	IAV
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction Influenza A virus subtype H5N1, also known as A(H5N1) or simply H5N1, is a subtype of the

influenza A virus which can cause illness in humans and many other animal species. A bird-adapted strain of H5N1, called HPAI A(H5N1) for highly pathogenic avian inf

Keywords

Nonstructural Protein 2; Nonstructural Protein 2 (H5N1)
