



Anti-Norovirus Capsid protein VP1 Polyclonal antibody (DPAB-DC4592)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Norovirus Capsid protein VP1
Immunogen	A synthetic peptide corresponding to N-terminus of Norovirus Capsid protein VP1 . The sequence is LTGQQNIIDPWIMNNFV
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Norovirus
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Norovirus is a genus of genetically diverse single-stranded RNA, non-enveloped viruses in the Caliciviridae family. The known viruses in the genus are all considered to be the variant strains
---------------------	---

of a single species called Norwalk virus. The viruses are trans

Keywords

Capsid protein VP1; Norovirus Capsid protein VP1
