



# Anti-NO-tyrosine polyclonal antibody (DPAB-DC4596)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |   |
|---------------------------|---|
| <b>Immunogen</b>          | A peptide (conjugated with KLH) corresponding to NO-tyrosine.       |
| <b>Isotype</b>            | IgG   |
| <b>Source/Host</b>        | Rabbit  |
| <b>Species Reactivity</b> | Human   |
| <b>Purification</b>       | Protein A purification  |
| <b>Conjugate</b>          | Unconjugated  |
| <b>Applications</b>       | WB,   |
| <b>Format</b>             | Liquid  |
| <b>Size</b>               | 100 µl  |
| <b>Buffer</b>             | In HEPES, 150 mM NaCl (50% glycerol, 0.01% BSA, 0.03% sodium azide) |
| <b>Preservative</b>       | 0.03% Sodium Azide  |
| <b>Storage</b>            | Store at -20°C. Aliquot to avoid repeated freezing and thawing.     |

## BACKGROUND

|                 |             |
|-----------------|-------------|
| <b>Keywords</b> | NO-tyrosine |
|-----------------|-------------|