



Anti-NPY (aa 1-20) polyclonal antibody (DPAB-DC4597)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with BSA) corresponding to amino acids 1-20 of cys (21)-Neuropeptide Y (NPY).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	IHC-Fr,
Format	Lyophilized
Size	50 µl
Buffer	Lyophilized from serum
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50-100 uL distilled water, store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Neuropeptide Y (NPY) is a 36-amino acid neuropeptide that acts as a neurotransmitter in the brain and in the autonomic nervous system of humans; slight variations of the peptide are found in many other animals. In the autonomic system it is produced mainly
---------------------	--

Keywords

NPY

GENE INFORMATION

Entrez Gene ID

[4852](#)

UniProt ID

[A4D158](#)
