



Anti-NTHL1 (full length) polyclonal antibody (DPAB-DC4598)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to full length NTHL1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Size	100 µl
Buffer	In antiserum (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at -20°C or -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	NTHL1
----------	-------

GENE INFORMATION

Entrez Gene ID [4913](#)

UniProt ID [E5KTI5](#)
