



Anti-Ovalbumin (native) polyclonal antibody (DPAB-DC4602)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified from chicken egg albumin, OVA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Chicken
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Size	100 µl
Buffer	In PBS (0.2% gelatin, 0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. Do not freeze.

BACKGROUND

Introduction	Ovalbumin (abbreviated OVA) is the main protein found in egg white, making up 60-65% of the total protein. Ovalbumin displays sequence and three-dimensional homology to the serpin superfamily, but unlike most serpins it is not a serine protease inhibitor.
Keywords	OVA; Ovalbumin