



Anti-Oxycodone polyclonal antibody (DPAB-DC4606)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Oxycodone conjugated with BSA.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Oxycodone is a semi-synthetic opioid synthesized from poppy-derived thebaine. It is a narcotic analgesic generally indicated for relief of moderate to severe pain. It was developed in 1916 in Germany as one of several new semi-synthetic opioids in an attempt to create a more potent and less addictive painkiller than morphine.
Keywords	Oxycodone