



Anti-OXT polyclonal antibody (DPAB-DC4609)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Oxytocin conjugated to bovine thyroglobuline.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Pig
Conjugate	Unconjugated
Applications	IHC-Fr,
Format	Lyophilized
Size	50 µl
Buffer	Lyophilized from serum
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50-100 uL distilled water, store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Oxytocin (Oxt) is a mammalian neurohypophysial hormone. Produced by the hypothalamus and stored and secreted by the posterior pituitary gland, oxytocin acts primarily as a neuromodulator in the brain.
Keywords	OXT; Oxytocin