



Anti-E. coli P2ox Polyclonal antibody (DPAB-DC4610)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	E. coli P2ox
Immunogen	Native purified Escherichia coli p2ox.
Isotype	IgG
Source/Host	Goat
Species Reactivity	E. coli
Conjugate	Unconjugated
Applications	WB (Recombinant protein), IP, ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	p2ox
----------	------