



Anti-PAF1 (full length) polyclonal antibody (DPAB-DC4611)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to full length drosophila Paf1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Flies, Fruit fly, Human
Conjugate	Unconjugated
Applications	WB, IF, IP,
Format	Liquid
Size	50 µg
Buffer	In PBS (0.05% sodium azide, 0.05% proclin 300)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Paf1
----------	------

GENE INFORMATION

Entrez Gene ID [54623](#)

UniProt ID [Q8N7H5](#)
