



Anti-Pancreastatin polyclonal antibody (DPAB-DC4612)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic Pancreastatin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr,
Format	Liquid
Size	250 µl
Buffer	In Tris-HCl buffer (1% BSA, 0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Chromogranin A is the precursor to several functional peptides including vasostatin, pancreastatin, catestatin and parastatin. These peptides negatively modulate the neuroendocrine function of the releasing cell (autocrine) or nearby cells (paracrine). Ot
Keywords	Pancreastatin

GENE INFORMATION

Entrez Gene ID [1113](#)

UniProt ID [P10645](#)
