



Anti-Pancreatic islets (native) polyclonal antibody (DPAB-DC4613)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross reactivity was determined by immunoblot (western) to be minimal against porcine kidney and porcine liver cell lysates.
Immunogen	Native purified pig pancreatic islets.
Source/Host	Dog
Species Reactivity	Pig
Conjugate	Unconjugated
Applications	WB, IP, ELISA,
Format	Liquid
Size	500 µg
Buffer	In 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The pancreas is a glandular organ in the digestive system and endocrine system of vertebrates. In humans it is located in the abdominal cavity behind the stomach. It is an endocrine gland producing several important hormones, including insulin, glucagon,
---------------------	--

