



Anti-Pantoea agglomerans polyclonal antibody (DPAB-DC4614)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross reactivity ratioEnterobacter cloacae : 1/5
Immunogen	Pantoea agglomerans total protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere.After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years.Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Pantoea agglomerans is a Gram-negative bacterium that belongs to the family
---------------------	---

Enterobacteriaceae.

Formerly called Enterobacter agglomerans, this bacterium is known to be an opportunistic pathogen in the immunocompromised, causing wound, blood, and urinary-t

Keywords

P. agglomerans; Pantoea agglomerans (Erwinia herbicola)
