



Anti-PCNA (aa 55-71) polyclonal antibody (DPAB-DC4617)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody would react with PCNA from the following alga sources: Prorocentrum minimum, Tetraselmis chui, Pleurochrysis carterae and Isochrysis galbana. Expect partial reactivity with PCNA from other alga sources as well as carrot, tobacco, peas, soybe
Immunogen	A synthetic peptide corresponding to amino acids 55-71 of Isochrysis galbana PCNA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB (Recombinant protein), IP, ELISA,
Format	Liquid
Size	100 µg
Buffer	In 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords PCNA

GENE INFORMATION

Entrez Gene ID [5111](#)

UniProt ID [P12004](#)
