



# Anti-PDS1 (internal region) polyclonal antibody (DPAB-DC4622)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Molecular weight: 42 KDa. Epitope: Internal region.
<b>Immunogen</b>	A synthetic peptide corresponding to internal region of human PDS1.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunoaffinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P,
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	100 µl
<b>Buffer</b>	In 10 mM PBS, pH 7.4 (1% BSA, 0.1% sodium azide)
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at 4°C. Do not freeze.

## BACKGROUND

