



Anti-PEPC (native) polyclonal antibody (DPAB-DC4623)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified PEPC from maize leaves.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Corn
Conjugate	Unconjugated
Applications	WB, IHC, ELISA,
Format	Liquid
Size	1 mg
Buffer	In 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	PEPC
-----------------	------