



Anti-Pepsin A polyclonal antibody (DPAB-DC4624)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Pepsin A polyclonal antibody detects endogenous levels of Pepsin A protein.
Immunogen	A synthetic peptide corresponding to Pepsin A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB (Cell lysate), IHC-P,
Format	Liquid
Concentration	1 mg/mL
Size	100 µl
Buffer	In PBS, pH 7.2 (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction

Pepsin A (EC 3.4.23.1, pepsin, lactated pepsin, pepsin fortior, fundus-pepsin, elixir lactate of pepsin, P I, lactated pepsin elixir, P II, pepsin R, pepsin D) is an enzyme. This enzyme catalyses the following chemical reaction.

Keywords

Pepsin A
