



Anti-PGF1-alpha polyclonal antibody (DPAB-DC4626)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This rabbit polyclonal antibody is specific to 6-Keto-Prostaglandin F1alpha (PGF1-alpha). Hapten Reactivity (%): 6-keto-PGF1alpha 100% PGE2 <1% PGF2alpha <1% PGD2 <1% PGE1 PGF1alpha <1% Thromboxane B2 <1%
Immunogen	6-Keto-Prostaglandin F1alpha (PGF1-alpha).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA, RIA,
Format	Liquid
Size	1 ml
Buffer	In PBS (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords PGF1-alpha