



Anti-PGF2-alpha polyclonal antibody (DPAB-DC4627)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity PGF2a 100%6-keto PGF1a 0.9%Thromboxane B2 16%13, 14-dihydro-15-keto PGF2a 0.06%13, 14-dihydro-6.15diketo PGF2a <0.01%PGD2 0.2%PGE1 0.1%PGE2 0.3%6-keto PGE1 0.03%PGF1a 21%Arachidonic acid <0.001%

Immunogen PGF2-alpha.

Source/Host Sheep

Species Reactivity N/A

Conjugate Unconjugated

Applications RIA,

Format Liquid

Size 100 µl

Buffer In antiserum

Preservative None

Storage Store at -80°C.Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords PGF2-alpha