



# Anti-PGP 9.5 polyclonal antibody (DPAB-DC4628)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to PGP 9.5. The sequence is ASSEDTLLKDAAKVCR
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB, IHC, IF,
Format	Liquid
Size	50 µl
Buffer	In whole serum (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	PGP 9.5
----------	---------