



# Anti-PHAP III (C-terminal) polyclonal antibody (DPAB-DC4629)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to C-terminus of human PHAP III.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.4 (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. Do not freeze.

## BACKGROUND

Keywords	PHAP III
----------	----------