



Anti-phenyl Acetic acid polyclonal antibody (DPAB-DC4632)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Phenyl acetic acid conjugated with carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	IHC, ICC, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	phenyl acetic acid
----------	--------------------