



# Anti-PHF2 (internal region) polyclonal antibody (DPAB-DC4635)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Reacts with residues 70-82 [KHGPGPTPDV KRVC] of the 120 KDa PHF protein.
<b>Immunogen</b>	A synthetic peptide corresponding to internal region of PHF2.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS (0.02% sodium azide)
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at 4°C for short term. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Keywords</b>	PHF2
-----------------	------

# GENE INFORMATION

Entrez Gene ID [5253](#)

---

UniProt ID [O75151](#)

---