



# Anti-Phytanic Acid polyclonal antibody (DPAB-DC4644)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Phytanic acid conjugated with carrier protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	N/A
<b>Purification</b>	Ammonium sulfate precipitation
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC,
<b>Format</b>	Lyophilized
<b>Size</b>	100 µl
<b>Buffer</b>	Lyophilized from PBS, pH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Phytanic acid (or 3,7,11,15-tetramethyl hexadecanoic acid) is a branched chain fatty acid that humans can obtain through the consumption of dairy products, ruminant animal fats, and certain fish. Western diets are estimated to provide 50-100 mg of phytani
---------------------	---

**Keywords**Phytanic acid

---