



Anti-Placental lactogen (native) polyclonal antibody (DPAB-DC4646)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified human Placental Lactogen (hPL).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunoaffinity purification
Conjugate	Unconjugated
Applications	IHC-P,
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	In 10 mM PBS, pH 7.4 (1% BSA, 0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Placental lactogen is a polypeptide placental hormone, part of the somatotropin family. Its structure and function is similar to that of growth hormone. It modifies the metabolic state of the
--------------	---

mother during pregnancy to facilitate the energy supply of the

Keywords	hPL; Placental Lactogen (hPL)
-----------------	-------------------------------

GENE INFORMATION

Entrez Gene ID	1442
-----------------------	----------------------

UniProt ID	A8K6C2
-------------------	------------------------
