



Anti-PP2C polyclonal antibody (DPAB-DC4648)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to Protein Phosphatase 2 C alpha/beta.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine, Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB (Tissue lysate), IP,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.08% sodium azide)
Preservative	0.08% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	PP2C
-----------------	------