



Anti-PRKD3 polyclonal antibody (DPAB-DC4653)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to human PRKD3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Size	100 µl
Buffer	In HEPES, 150 mM NaCl (50% glycerol, 0.01% BSA, 0.03% sodium azide)
Preservative	0.03% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Serine/threonine-protein kinase D3 or protein kinase C $\delta 3$ (PKC- $\delta 3$) is an enzyme that in humans is encoded by the PRKD3 gene.
Keywords	PRKD3

GENE INFORMATION

Entrez Gene ID [23683](#)

UniProt ID [O94806](#)
