



# Anti-Procollagen Type III (C-terminal) polyclonal antibody (DPAB-DC4654)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	This antibody reacts with human and bovine procollagen type III (PIIIP) 100%; human and bovine collagen type III very low at RIA dilution <1:100. human and bovine procollagen type I <0.1% in RIA at 1:200 dilution.
<b>Immunogen</b>	Procollagen type III deprived of the C-terminus PIIICP and N-terminus PIINP.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Bovine, Human, Pig
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P, IF, ELISA, RIA,
<b>Format</b>	Lyophilized
<b>Size</b>	100 µl
<b>Buffer</b>	No additive
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C. Aliquot after reconstitution to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Procollagen peptidase (EC 3.4.24.14, procollagen N-terminal peptidase, procollagen aminopeptidase, aminoprocollagen peptidase, aminoterminal procollagen peptidase,
---------------------	---

procollagen aminoterminal protease, procollagen N-terminal proteinase, type I/II procollage

---

**Keywords**

Procollagen Type III

---