



Anti-Progesterone polyclonal antibody (DPAB-DC4655)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Progesterone (4)-BTG.
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	EIA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Progesterone (pregn-4-ene-3,20-dione; abbreviated as P4) is an endogenous steroid hormone involved in the menstrual cycle, pregnancy, and embryogenesis of humans and other species. It belongs to a group of steroid hormones called the progestogens, and is
Keywords	P4; Progesterone