



# Anti-Pro-KLK5 polyclonal antibody (DPAB-DC4657)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	This antibody recognize both native and denatured KLK5.
Immunogen	Recombinant protein corresponding to pro-kallikrein 5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB,
Format	Lyophilized
Size	200 µl
Buffer	Lyophilized from serum
Preservative	None
Storage	Store at -20°C on dry atmosphere. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	Pro-KLK5
----------	----------

## GENE INFORMATION

Entrez Gene ID [25818](#)

---

UniProt ID [Q9Y337](#)

---