



# Anti-Propoxyphene polyclonal antibody (DPAB-DC4659)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Propoxyphene conjugated with BSA.
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	N/A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	EIA,
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Dextropropoxyphene is an analgesic in the opioid category, patented (1955) and manufactured by Eli Lilly and Company. It is intended to treat mild pain and also has antitussive (cough suppressant) and local anaesthetic effects. The drug has been taken off
<b>Keywords</b>	Propoxyphene