



# Anti-Protein C (native) polyclonal antibody (DPAB-DC4661)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Native purified Protein C from bovine plasma.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Bovine
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, EIA,
<b>Format</b>	Liquid
<b>Size</b>	150 µg
<b>Buffer</b>	In PBS, pH 7.5 (0.02% sodium azide, 50% glycerol (v/v))
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Protein C, also known as autoprothrombin IIA and blood coagulation factor XIV,:6822 is a zymogenic (inactive) protein, the activated form of which plays an important role in regulating blood clotting, inflammation, cell death, and maintaining the permeabi
---------------------	---

