



Anti-Protein Kinase Consensus Sequence polyclonal antibody (DPAB-DC4663)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It is excellent for Western blot detection of protein kinases containing this sequence.
Immunogen	A synthetic peptide corresponding to Protein Kinase Consensus Sequence. The sequence is GTPEYLAPE
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	1 ml
Buffer	In PBS
Preservative	None
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Protein Kinase Consensus Sequence
-----------------	-----------------------------------