



Anti-Protein L polyclonal antibody (DPAB-DC4664)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to Protein L.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB (Recombinant protein),
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (30% glycerol, 0.5 mg/mL BSA, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Protein L was first isolated from the surface of bacterial species <i>Peptostreptococcus magnus</i> and was found to bind immunoglobulins through L chain interaction, from which the name was suggested. It consists of 719 amino acid residues. The molecular weight
Keywords	Protein L