



# Anti-PRPH polyclonal antibody (DPAB-DC4666)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |   |
|---------------------------|---|
| <b>Immunogen</b>          | Recombinant protein corresponding to PRPH.                      |
| <b>Isotype</b>            | IgG   |
| <b>Source/Host</b>        | Rabbit  |
| <b>Species Reactivity</b> | Mammals   |
| <b>Conjugate</b>          | Unconjugated  |
| <b>Applications</b>       | WB, ICC, IF,  |
| <b>Format</b>             | Liquid  |
| <b>Size</b>               | 100 µl  |
| <b>Buffer</b>             | In PBS (50% glycerol, 0.01% thimerosal)                         |
| <b>Preservative</b>       | None  |
| <b>Storage</b>            | Store at -20°C. Aliquot to avoid repeated freezing and thawing. |

## BACKGROUND

|                 |      |
|-----------------|------|
| <b>Keywords</b> | PRPH |
|-----------------|------|