



Anti-*P. aeruginosa* Polyclonal antibody (DPAB-DC4668)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts predominately against <i>P. aerug.</i> with some cross reactivity against other <i>Pseudomonas</i> species (<i>P. maltophilia</i> , <i>P. fluorescens</i> and <i>P. gladioli</i>).
Target	<i>P. aeruginosa</i>
Immunogen	Native purified from outer membrane protein extract of <i>Psuedomonas aeruginosa</i> .
Source/Host	Pig
Species Reactivity	<i>P. aeruginosa</i>
Conjugate	Unconjugated
Applications	sELISA, Agg, Detection Antibody,
Format	Liquid
Size	1 ml
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	<i>Pseudomonas aeruginosa</i> is a common bacterium that can cause disease in animals, including humans. It is citrate, catalase, and oxidase positive. It is found in soil, water, skin flora, and
---------------------	---

most man-made environments throughout the world. It thrives not on

Keywords

P. aeruginosa; *Pseudomonas aeruginosa*
