



Anti-*P. aeruginosa* Polyclonal antibody (DPAB-DC4671)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	<i>P. aeruginosa</i>
Immunogen	<i>Pseudomonas aeruginosa</i> total protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>P. aeruginosa</i>
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	WB,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction *Pseudomonas aeruginosa* is a common bacterium that can cause disease in animals, including

humans. It is citrate, catalase, and oxidase positive. It is found in soil, water, skin flora, and most man-made environments throughout the world. It thrives not on

Keywords	P. aeruginosa; Pseudomonas aeruginosa
-----------------	---------------------------------------
