



# Anti-*P. aeruginosa* Tail Fiber Polyclonal antibody (DPAB-DC4674)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	<i>P. aeruginosa</i> Tail Fiber
<b>Immunogen</b>	A synthetic peptide corresponding to N-terminus of <i>Pseudomonas aeruginosa</i> Tail Fiber.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	<i>P. aeruginosa</i>
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
<b>Preservative</b>	0.01% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	<i>Pseudomonas aeruginosa</i> is a common bacterium that can cause disease in animals, including humans. It is citrate, catalase, and oxidase positive. It is found in soil, water, skin flora, and most man-made environments throughout the world. It thrives not on
---------------------	--

**Keywords**

Tail Fiber; *Pseudomonas aeruginosa* Tail Fiber

---